

FERROSOUND

## PL125-D30 PERSONAL AUDIO LOOP

SOUND INNOVATIONS FOR HEARING IMPAIRED PEOPLE

### About the Ferrosound PL125-D30

Personal Induction Loops may be used wherever a normally hearing person would use headphones.

This new version of the PL125 is designed for use with TV and PC headphone outputs, where current may be limited to prevent headphone volumes from being too loud.

The PL125-D is fully stereo compatible, in which both left and right channels are mixed within the loop. This 'D' version will also handle higher currents, making it more suitable for use with powerful conversation aids, such as the Ferrosound System 121.

The PL125-D is available in two versions. You may choose either a loop with an integral lead or a loop which plugs into an external lead incorporating a volume control.

Both versions are fitted with a 3.5mm stereo plug and supplied with an adaptor for use with 6.25mm sockets.



*Volume control with 0.9m lead*



### Specifications

Impedance @ 1kHz	30 + 30 ohms (nominal)
Frequency response *	100 Hz –8kHz $\pm$ 3dB
Loop length	900mm
Lead length	0.9 m
Option 1	Case fitted with 3.5mm stereo plug for use with 0.9m lead & volume control
Option 2	Integral lead fitted with 3.5mm plug.

\* Frequency response and sensitivity of induction loops will vary according to the electrical characteristics of the audio source used.